




<p>CVS Centrifuge Tube: <b>Model # 45170</b>, 12.5 mL. <b>Finger Tube</b>          Made to API Standards, and graduated as follows:</p> <ul style="list-style-type: none"> <li>• 0 - 3% in 0.2% Graduations</li> <li>• 3 - 10 % in 0.5% Graduations</li> <li>• 10 - 50% in 1.0% Graduations</li> <li>• 100% equals 12.5 mL</li> </ul> <p>Finger Tube Shield - Part Number CVS18371</p>	
<p>CVS Centrifuge Tube: <b>Model # 45244 - 100</b>, 100mL <b>Pear Shaped</b>          Made to API Standards, and graduated as follows:</p> <ul style="list-style-type: none"> <li>• 0 - 1.5 mL in 0.1 mL Graduations</li> <li>• 1.5 - 5 mL in 0.5 mL Graduations</li> <li>• 5 - 10 mL in 1.0 mL Graduations</li> <li>• 15, 20, 25, 50, 100 mL Graduations</li> </ul> <p>CVS Centrifuge Tube: <b>Model # 45244-200</b>, 200% <b>Pear Shaped</b>          Made to API Standards and graduated as follows:</p> <ul style="list-style-type: none"> <li>• 0 - 3% in 0.1% Graduations</li> <li>• 3 - 5 % in 0.5% Graduations</li> <li>• 5 - 10 % in 1.0% Graduations</li> <li>• 10 - 100% in 10% Graduations</li> <li>• 100 - 200% in 20% Graduations</li> </ul> <p>Pear Shaped Shield - Part Number CVS18080-R</p>	
<p>CVS Centrifuge Tube: <b>Model # 45243- 100</b>, 100mL <b>Short Cone</b>          Made to API Standards, and graduated as follows:</p> <ul style="list-style-type: none"> <li>• 0 - 0.5 mL in 0.05 mL Graduations</li> <li>• 0.5 - 2 mL in 0.1 mL Graduations</li> <li>• 2 - 3 mL in 0.2 mL Graduations</li> <li>• 3 - 5 ml in 0.5 mL Graduations</li> <li>• 5 - 10 mL in 1.0mL Graduations</li> <li>• 10 - 25 mL in 5mL Graduations</li> <li>• 50, 75, 100 mL Graduations</li> </ul> <p>CVS Centrifuge Tube: <b>Model # 45243-200</b>, 200% <b>Short Cone</b>          Made to API Standards and graduated as follows:</p> <ul style="list-style-type: none"> <li>• 0 - 1% in 0.1% Graduations</li> <li>• 1 - 4 % in 0.2% Graduations</li> <li>• 4 - 6 % in 0.4% Graduations</li> <li>• 6 - 10% in 1.0% Graduations</li> <li>• 10 - 20% in 2.0% Graduations</li> <li>• 20 - 50% in 10.0% Graduations</li> <li>• 100%, 150%, 200% Graduations</li> </ul> <p>Short Cone Shield - Part Number CVS18081R</p>	
<p>CVS Centrifuge Tube: <b>Model # 45241- 100</b>, 100mL <b>Long Cone</b>          Made to API Standards          (Same graduations as the CVS Short Cone 100 ml Tubes)          CVS Centrifuge Tube: <b>Model # 45241-FG</b>, 100mL <b>Long Cone</b>          Fully Graduated</p>	

# CVS

---

## Controls Ltd.

### **CVS Controls Ltd. Centrifuge Tube Holder**

The CVS Centrifuge Tube Holder allows you to easily manage you samples, saving time and money due to damaged tubes or lost samples.

Designed to accept both Shortcone, or Pear Shaped tubes.

Standard with 4 slots.

Can easily be wall mounted or placed on workbench. Rubber feet keep the CVS Tube Holder from slipping on flat surfaces.

Durable and long lasting hammer tone powder coating.

**Part Number:**

CVS S10D-RACK—Shortcone Tubes

CVS S10B-RACK—Pear Shaped Tubes



For more information contact a CVS Controls Representative or visit us online at [www.cvs-controls.com](http://www.cvs-controls.com)